

LABORHITS

Reflectance Spectrophotometer TR 520



Features

- Dual apertures
- Measurement of fluorescent materials with UV on / off
- Extended range of illuminants

Specifications

Light source: combined LED light, UV light

Illuminant: D65, A, C, D50, D55, D75, F1, F2, F3, F4, F5, F6, F7, F8, F9, F10, F11, F12

Blende: 10/8 + 5/4 mm

Repeatability: MAV/SCI: $\Delta E \leq 0,03$

Inter instrument agreement:

MAV/SCI: $\Delta E \leq 0,15$

Dimensions (W x D x H): 184 x 105 x 77 mm

Weight: 600 g

Scope of supply:

Device, OnShade PC software, black and white calibration board (white tile and black interface), dust cover, USB cable, external DC power supply, USB stick with user manual in the operating languages, carrying case



Type	PK	Cat. No.
TR 520	1	4.672 487