

# LABORHITS

## Roller Mixers LLG-uniRoller 6 easy and pro

For gentle rolling of samples

Max speed 70 rpm

Easy to clean



**LLG-uniRoller 6 easy**

- Analog control
- Easy to use



**LLG-uniRoller 6 pro**

- LCD display
- Continuous or timer operation



Analog or digital tube rollers with 6 rolls for the efficient mixing of biological samples, viscous substances and liquid/solid suspensions. Used in a variety of applications, such as the prevention of blood coagulation, immune precipitation.

**Specifications**

	LLG-uniROLLER 6	
	easy	pro
Type of movement:	rocking and rolling	
Amplitude:	24 mm	
Speed:	0 ... 70 rpm	
Timer:	-	1 ... 1199 min
Load max.:	4 kg	
Ambient temperature:	5 ... 40 °C	
Max. humidity:	80 % rel. F.	
Dimensions (W x D x H):	260 x 450 x 120 mm	
Weight:	4.5 kg	5.1 kg
Power consumption:	25 W	30 W
Power supply:	100 ... 240 V, 50/60 Hz	
IP code:	IP 21	

Type	Plug type	PK	Cat. No.
LLG-uniROLLER 6 easy	EU	1	<b>6.263 645</b>
LLG-uniROLLER 6 easy	UK	1	<b>6.263 646</b>
LLG-uniROLLER 6 pro	EU	1	<b>6.263 647</b>
LLG-uniROLLER 6 pro	UK	1	<b>6.263 648</b>