

LABORHITS

HPLC Reservoir Bottles DURAN®

For longer, uninterrupted analyses in HPLC

DURAN® borosilicate 3.3 glass

With conical base

NEW



Features

- GL 45 Thread
- ISO 3585
- Type I neutral glass, by USP / EP
- With Retrace Code and certificate for manufacturing lot traceability
- Available in four sizes: 1, 2, 5, or 10 Litres
- 5 and 10 litre bottles can accommodate magnetic stirring bars
- Durable white enamel volume graduations
- Autoclavable at 121 °C or 134 °C

The DURAN® reservoir bottles have been designed for use in high-performance liquid chromatography (HPLC) and allows longer, uninterrupted analyses. A central cavity in the bottom of the bottle allows mobile phases and solvents to be dispensed without tilting the bottle. Side vent hole in the base for ultrasonic bath degassing, or water drainage after cleaning. Screw caps HPLC (PTFE, GL 45, 3 ports) available as accessories (Art.-No.: 7.920 475).

Nominal capacity ml	Diam. mm	Height mm	PK	Cat. No.
1,000	110	295	1	4.665 543
2,000	145	309	1	4.665 544
5,000	190	386	1	4.665 545
10,000	235	481	1	4.665 546