

LABORHITS

Midi Centrifuge G-L

For all types of separation processes

High-performance motor



Features

- With metal cover
- Starts only when the cover is firmly closed
- Direct and easy setting of speed, rotational speed and programmes via membrane keypad

Due to the material separation at medium rotational speeds, this midi centrifuge is ideal for all types of separation processes, such as in molecular biology. It can be used with standard 1.5 / 2 ml tubes as well as with PCR strips. The cover opens automatically at the press of a button.

Scope of delivery: Midi centrifuge incl. metal rotor and cover for 12 vessels, with a capacity of 1.5/2 ml

Specifications

Max. RCF: 16,500 x g
Max. speed: 15,700 min⁻¹
Max. capacity: 12 x 1.5 / 2.0 ml
Power consumption: 65 W
Dimensions (W x D x H): 200 x 280 x 140 mm
Weight: 3.9 kg
Power supply: 220 ... 230 V, 50 ... 60 Hz
Protection class: IP 20

Type	PK	Cat. No.
G-L	1	4.664 755