



BRAND

Reagecon

Delivering the Correct Result



Save now
up to
38%*

accu-jet® pro

Don't miss our special price for Pipetting aid accu-jet® pro only from September 1 until December 31, 2020!

BRAND. For lab. For life.®

Pipetting aid accu-jet® pro

All with one hand: select the delivery mode (gravity-delivery / blow-out) and adjust the motor speed range with your thumb; use variable button pressure for finest control of filling and delivery speed. The ergonomic handgrip, perfectly balanced design, weight only 190 g - for fatigue-free pipetting even in prolonged operations.

The accu-jet® pro is supplied with two battery compartment covers. The second cover with ridges enables you to place the instrument inverted on the bench, to avoid contamination of the pipette adapter.

Your benefits:

- + Powerful: Fills a 50 ml pipette in less than 10 seconds.
- + Sensitive: Precise selection of pipetting speed.
- + Long-lasting: up to 8 hours of continuous pipetting without recharging

Promo accu-jet® pro

Promotion runs from September 1 until December 31, 2020

Save up to 38%*!

Simply enter the desired quantity in the table and send the completed form to your distributor.



Yes, I would like to participate in the promotion and order at a special price*:

accu-jet® pro Pipetting aid

with nickel-metal hydride battery, 2 battery compartment covers, wall support, AC adapter (100-240 V; 50/60 Hz) and 2 spare membrane filters 0.2 µm, sterile

Cat. No.	Description	Quantity	Special price €*	Total price €*
4000031	dark blue, AC adapter for Ireland/UK 100-240 V; 50/60 Hz		249,00	
4000032	magenta, AC adapter for Ireland/UK 100-240 V;50/60 Hz		249,00	
4000033	green, AC adapter Ireland/UK 100-240 V; 50/60 Hz		249,00	
4000034	royal blue, AC adapter Ireland/UK 100-240 V;50/60Hz		249,00	

Total:

Date

Signature

Address:

.....
First name Last name

.....
Institute / Company

.....
Street Zip code and City

.....
E-Mail Tel./Fax

Distributor:

Reagecon
Delivering the Correct Result

* All prices are recommended sales prices without VAT. Savings in percent compared with BRAND price list 2020.

Your order will be delivered through a distributor.

BRAND®, BRAND. For lab. For life.®, accu-jet® S as well as the figurative marks depicted here and the BRAND figurative mark are registered trademarks or trademarks of BRAND GMBH + CO KG, Germany. All other trademarks mentioned or depicted here are the property of the respective owners. Subject to technical modification without notice. Errors excepted.

Reagecon
Delivering the Correct Result

Reagecon Diagnostics Ltd.
Shannon Free Zone
Shannon, Co. Clare
Ireland

Tel: +353 61 472622
Fax: +353 61 472642
Mail: sales@reagecon.ie
Web: www.reagecon.com

BRAND