

# LABORHITS

## Digital thermometer testo 720

**Ideal for Use in Laboratories**

**Versatile and Precise**



### testo 720

- Connection for Pt100 sensor
- Large illuminated display
- Continuous display of min./max. values
- Audible warning signal
- Hold-function

**Scope of supply:** testo 720, calibration protocol, batteries

**Please order sensor seperately.**

### TopSafe

The solid, waterproof protective cover protects the measuring device from moisture, impact and dirt.

### Specifications

**Measuring range:** -100 to +800 °C

**Operating temperature:** -20 to +50 °C

**Housing material:** ABS

**Accuracy:** System measurements with immersion / penetration probes  
±0.2 % f.s.d. (+200 ... +800 °C)  
±0.2 °C (remaining range)

**Resolution** 0.1 °C

**Battery type:** 9 V block battery

**Dimensions (L x W x H):** 182 x 64 x 40 mm

**Weight:** 171 g

Type	PK	Cat. No.
testo 720	1	<b>9.851 323</b>
Laboratory probe Pt100, glass-coated, glass tube DURAN®	1	<b>9.851 718</b>
TopSafe protective cover against impact and dust	1	<b>9.718 020</b>