

LABORHITS

Laboratory Bottles Protect DURAN®

Resistant , Transparent Plastic Coating

DURAN® Borosilicate Glass 3.3

With Retrace Code

Highly transparent. The coating provides scratch, leak and splinter protection and is ideally suited to both the transport and storage of hazardous media or valuable samples. Service temperature limit of the PU plastic coating: -30 °C to 135 °C. With screw cap and pouring ring.

Features

- Graduated
- UV protection
- DIN thread
- DIN ISO 4796-1
- Suitable for microwaving
- USP compliant
- Autoclavable



Nominal capacity ml	Outer diam. mm	Height mm	Thread GL	PK	Cat. No.
100	56	100	45	1	4.652 768
150	62	110	45	1	4.652 769
250	70	138	45	1	4.652 770
500	86	176	45	1	9.071 990
1000	101	225	45	1	9.071 994
2000	136	260	45	1	9.072 006
3500	160	300	45	1	4.652 772
5000	182	330	45	1	4.652 773