

# LABORHITS

## Vortex Mixers ZX4/WIZARD/TX4 with IR Sensor

ZX4



WIZARD



TX4



### Resistant and Versatile

Up to 3000 rpm

3 Years Warranty

#### Features

- For IR sensor and continuous operation
- Electronic speed control
- Adjustable shaking speeds up to 3000 min<sup>-1</sup>
- Very good chemical resistance
- Broad range of accessories available

Vortexers with patented IR sensor for sample recognition: the shaking process starts automatically when the sample tube approaches the mixing cup head. This enables fatigue-free working, as the sample doesn't have to be pressed on the device.

**ZX4:** Analog vortex mixer

**WIZARD:** Compact, analog vortex mixer

**TX4:** Digital vortex mixer with timer and LCD display

**Scope of delivery:** Vortex mixer with mixing head

#### Specifications

**Type of movement:** orbital

**Orbital diameter:** 4.5 mm

**Speed range:** 0 ... 3000 min<sup>-1</sup>

**Input power:** 15 W

**IP code:** IP 42

**Power supply:** 100 ... 240 V, 50 ... 60 Hz

**Warranty:** 3 years

Type	Weight kg	Plug type	PK	Cat. No.
ZX4	2.7	EU	1	<b>7.971 461</b>
ZX4 with foam stands	2.7	EU	1	<b>4.666 714</b>
WIZARD	2.2	EU	1	<b>4.666 536</b>
TX4	2.7	EU	1	<b>4.666 537</b>