

LABORHITS

Heating Ovens Series FD



Excellent Temporal and Spatial
Temperature Accuracy

Temperature Range up to +300 °C



Features

- Temperature range: room temperature plus 10 °C to 300 °C
- Controller with LCD display
- Electromechanical control of the exhaust air flap
- Class 2 independent temperature safety device (DIN 12880) with visual temperature alarm
- USB port for recording data
- 2 grids included
- 230 V, 50/60 Hz

The drying chambers with forced convection are renowned due to their new design, convenient operation and efficiency. The latest APT.line™ technology offers outstanding temperature accuracy.

Type	Nominal capacity l	Rated capacity kW	Internal dimensions (W x D x H) mm	External dimensions (W x D x H) mm	Max. grids	Weight kg	PK	Cat. No.
FD 56	60	1.10	400 x 345 x 440	560 x 565 x 625	4	39	1	6.280 069
FD 115	116	1.30	550 x 385 x 550	710 x 605 x 735	5	50	1	4.658 060