

LABORHITS

OBE • OBE-S • OBS

KERN OBE-11



KERN OBE 112

KERN OBE-S



KERN OBE 114T241

KERN OBS-10



KERN OBS 106

KERN OBE-11

- Strong and continuously dimmable 3 W LED
- Focusable 1.25 ABBE condenser
- 4-fold lens revolving unit
- Siedentopf 30° inclined, 360° rotatable

KERN OBE 114T241

- 4-fold lens revolving unit
- Siedentopf 30° inclined, 360° rotatable
- Dioptr compensation one-sided
- Set of microscope and tablet

KERN OBS-10

- Strong and continuously dimmable 0.5 W LED
- Focusable 1.25 ABBE condenser
- 4-fold lens revolving unit
- Siedentopf 30° inclined, 360° rotatable

Model	Standard configuration					Cat. No.
	Tube	Eyepiece	Lens quality	Lens	Illumination	
KERN OBE 112	Binocular	HWF 10×/Ø 18 mm	Achromatic	4× / 10× / 40× / 100×	3 W LED	7.940 489
OBE 114	Trinocular	HWF 10×/Ø 18 mm	Achromatic	4× / 10× / 40× / 100×	3 W LED	6.285 896
OBE 114T241	Binocular	HWF 10×/Ø 18 mm	Achromatic	4× / 10× / 40× / 100×	3 W LED	6.318 003
OBS 105	Monocular	WF 10×/Ø 18 mm	Achromatic	4× / 10× / 40×	0,5 W LED	4.681 238
OBS 106	Binocular	WF 10×/Ø 18 mm	Achromatic	4× / 10× / 40×	0,5 W LED	6.285 895

LABORHITS

OBL • OZL

KERN OBL-12•13



KERN OBL 137

KERN OZL-46



KERN OZL 466

KERN OBL-12•13

- Strong and continuously dimmable 3 W LED
- Focusable 1.25 ABBE condenser
- 4-fold lens revolving unit
- Dioptre compensation one-sided
- Siedentopf 30° inclined, 360° rotatable

KERN OZL-46

- Strong, continuously dimmable, integrated LED ring illumination
- Brightness adjustable (separate)
- Magnification ratio: 6,4:1
- Dioptre compensation both sides
- Lens barrel 45° tilt
- The zoom objective offers you continuous magnification from 7× – 45×

Model	Standard configuration				Cat. No.
	Tube	Eyepiece	Lens	Illumination	
KERN OBL 127	Binocular	HWF 10×/Ø 20 mm	4× / 10× / 40× / 100×	3 W LED	4.657 659
OBL 137	Trinocular	HWF 10×/Ø 20 mm	4× / 10× / 40× / 100×	3 W LED	4.657 661
OZL 463	Binocular	HWF 10×/Ø 20 mm	0,7× - 4,5×	3 W LED	6.286 167
OZL 466	Binocular	HWF 10×/Ø 20 mm	0,7× - 4,5×	3 W LED	4.679 951