

LABORHITS

Compact balances Compass™ CR



Slim Stackable Design

Large Weighing Platform

Easy to Use



Features

- Battery life up to 300 hours in continuous operation
- Backlit LCD display, digit height 15 mm
- Easy and convenient 2-button operation
- ABS housing and weighing platform with non-slip feet

Economic balance series for basic weighing. With a large weighing platform for a variety of mobile weighing applications in laboratories, industry and education.

Scope of supply: Balance incl. 3 AA batteries included.

Please order power adapter separately.

Specifications

Dimensions (W x D x H): 140 x 205 x 41 mm

Dimensions plate (W x D): 125 x 132 mm

Weight: 400 g

Type	Capacity g	Readability g	PK	Cat. No.
CR221	220	0,1	1	6.275 420
CR621	620	0,1	1	6.275 421
CR2200	2200	1,0	1	6.275 422
CR5200	5200	1,0	1	6.275 423