

LABORHITS

Precision Balances EMB



Entry level laboratory balance with tremendous weighing performance

AA
QUALITY

Features

- Simple and convenient 2-key operation
- Tare function facilitates formulation work
- Particularly flat design
- Ready for use: Batteries included
- Ring-shaped draught shield standard, only for models with weighing plate size Ø 82 mm
- Hook for underfloor weighing standard
- Large LCD display, digit height 15 mm
- AUTO-OFF function after 3 min

The Precision balances EMB are entry-level devices in the field of laboratory balances, and they cover basic applications. They feature simple operation and the very best price/performance ratio. In addition they can be stacked and can also be operated by battery, as an alternative to mains power, which means that they have the flexibility to be used in different locations.

Specifications

Dimensions (W x D x H): 170 x 240 x 39 mm
EMB 100-3/EMB1000-2: 170 x 240 x 54 mm
Net weight: approx. 750 g

Type	Weighing plate Ø mm	Weighing capacity g	Readability g	Reproducibility g	Linearity ± g	PK	Cat. No.
EMB 100-3	82	100	0.001	0.001	0.003	1	6.238 426
EMB 200-2	105	200	0.01	0.01	0.02	1	9.901 320
EMB 600-2	105	600	0.01	0.01	0.03	1	9.901 321
EMB 1000-2	150	1000	0.01	0.01	0.05	1	6.235 728
EMB 2000-2	150	2000	0.01	0.01	0.05	1	7.670 022
EMB 500-1	150	500	0.1	0.1	0.2	1	9.901 322
EMB 1200-1	150	1200	0.1	0.1	0.3	1	9.901 323
EMB 5.2K1	150	5200	1.0	1.0	3.0	1	9.901 324
EMB 5.2K5	150	5200	5.0	5.0	10.0	1	9.901 372
EMB 2200-0	150	2200	1.0	1.0	2.0	1	9.901 371