

LABORHITS

Battery-powered samplers for liquids



Thanks to a suction depth of 5 meters, the battery-powered Vampire sampler easily pumps samples, even from sampling points at lower depths, for example from barrels, tanks, wells, water traps or wastewater and sewage treatment plants.

- Battery-powered hose pump for sampling liquids from depth
- Easy cleaning and exchangeable hoses prevent impurities
- Suction depth up to 5 metres
- Consisting of a drive unit, pump head, pump hoses (30 cm, Marprene/silicone), suction hose (5 m, PE), hose connector (PP), sample bottle (250 ml, LDPE), bottle basket, two rechargeable batteries and battery charger

Moreover the portable sampler is easy to handle and carry.



Type	PK	Cat. No.
Vampire sampler	1	9.303 851

LABORHITS

Battery-powered samplers for powders and granules



Quality assurance starts with the inspection of incoming goods! An all-layer sample of bulk goods is no longer a problem with PowderProof®, the sampler for powders and granules.

- Many uses
- Stainless steel AISI 316L (1.4404)/AISI 301 (1.4310)
- Incl. motor, two rechargeable batteries, fast charger, lance, spiral screw, bags, clip for holding the bag and cleaning brush
- Additional spiral screws available as accessory



Model	Lance	Insertion depth mm	PK	Cat. No.
for powders up to 2 mm dia.	horizontal	300	1	4.687 799
for powders up to 2 mm dia.	vertical	600	1	4.687 800
for powders up to 2 mm dia.	vertical	900	1	6.250 303
for powders up to 2 mm dia.	vertical	1200	1	6.262 872
for granules from 2 up to 4 mm dia.	horizontal	300	1	6.290 980
for granules from 2 up to 4 mm dia.	vertical	600	1	4.687 801
for granules from 2 up to 4 mm dia.	vertical	900	1	6.223 346
for granules from 2 up to 4 mm dia.	vertical	1200	1	4.687 802